

PIPE-T, PIPE-T8, PIPE-T12, PIPE-TF

FEATURES (SINGLE CHANNEL & MULTI-CHANNEL)	
ITEM	DESCRIPTION
1	Pipettes cover a volume range from 0.1µL to 10mL.
2	Ergonomic design provides excellent operating experience.
3	Large display window allows for easy volume identification.
4	Easy calibration and maintenance.
5	Online calibration is available.

FEATURES (MULTI-CHANNEL)	
ITEM	DESCRIPTION
1	PIPE-T8 and PIPE-T12 channel pipettes are appropriate for 96 well plates
2	Dispensing head rotates for effortless pipetting convenience.
3	Individual piston and tip cone assemblies allowing easy repair and maintenance.
4	Compound material-made tip cone secures high sealing performance.
5	Compatible with most universal tip brands.
6	Online calibration is available.



PIPE - T



PIPE - T8



PIPE - T12



PIPE - TF

All rights reserved by iSiinitzan. Reproduction and distribution of any part of this document is prohibited.



PIPE-T, PIPE-T8, PIPE-T12, PIPE-TF

TECHNICAL SPECIFICATION PIPE-T SINGLE-CHANNEL ADJUSTABLE VOLUME PIPETTES						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
0.1-2.5µL	0.05µL	2.5µL	2.50%	0.0625	2.00%	0.05
		1.25µL	3.00%	0.0375	3.00%	0.0375
		0.25µL	12.00%	0.03	6.00%	0.015
0.5-10µL	0.1µL	10µL	1.00%	0.1	0.80%	0.08
		5µL	1.50%	0.075	1.50%	0.075
		1µL	2.50%	0.025	1.50%	0.015
2-20µL	0.5µL	20µL	0.90%	0.18	0.40%	0.08
		10µL	1.20%	0.12	1.00%	0.1
		2µL	3.00%	0.06	2.00%	0.04
5-50µL	0.5µL	50µL	0.60%	0.3	0.30%	0.15
		25µL	0.90%	0.225	0.60%	0.15
		5µL	2.00%	0.1	2.00%	0.1
10-100µL	1µL	100µL	0.80%	0.8	0.15%	0.15
		50µL	1.00%	0.5	0.40%	0.2
		10µL	3.00%	0.3	1.50%	0.15
20-200µL	1µL	200µL	0.60%	1.2	0.15%	0.3
		100µL	0.80%	0.8	0.30%	0.3
		20µL	3.00%	0.6	1.00%	0.2
50-200µL	1µL	200µL	0.60%	1.2	0.15%	0.3
		100µL	0.80%	0.8	0.30%	0.3
		50µL	1.00%	0.5	0.40%	0.2
100-1000µL	5µL	1000µL	0.60%	6	0.20%	2
		500µL	0.70%	3.5	0.25%	1.25
		100µL	2.00%	2	0.70%	0.7
200-1000µL	5µL	1000µL	0.60%	6	0.20%	2
		500µL	0.70%	3.5	0.25%	1.25
		200µL	0.90%	1.8	0.30%	0.6
1000-5000µL	50µL	5000µL	0.50%	25	0.15%	7.5
		2500µL	0.60%	15	0.30%	7.5
		1000µL	0.70%	7	0.30%	3
2-10mL	0.1mL	10mL	0.60%	60	0.20%	20
		5mL	1.20%	60	0.30%	15
		2mL	3.00%	60	0.60%	12

TECHNICAL SPECIFICATION PIPE-T8-CHANNEL ADJUSTABLE VOLUME PIPETTES						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
0.5-10µL	0.1µL	10µL	1.50%	0.15	1.50%	0.15
		5µL	2.50%	0.125	2.50%	0.125
		1µL	4.00%	0.04	4.00%	0.04
5-50µL	0.5µL	50µL	1.00%	0.5	0.50%	0.25
		25µL	1.50%	0.375	1.00%	0.25
		5µL	3.00%	0.15	2.00%	0.1
50-300µL	5µL	300µL	0.70%	2.1	0.25%	0.75
		150µL	1.00%	1.5	0.50%	0.75
		50µL	1.50%	0.75	0.80%	0.4

TECHNICAL SPECIFICATION PIPE-T12-CHANNEL ADJUSTABLE VOLUME PIPETTES						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
0.5-10µL	0.1µL	10µL	1.50%	0.15	1.50%	0.15
		5µL	2.50%	0.125	2.50%	0.125
		1µL	4.00%	0.04	4.00%	0.04
5-50µL	0.5µL	50µL	1.00%	0.5	0.50%	0.25
		25µL	1.50%	0.375	1.00%	0.25
		5µL	3.00%	0.15	2.00%	0.1
50-300µL	5µL	300µL	0.70%	2.1	0.25%	0.75
		150µL	1.00%	1.5	0.50%	0.75
		50µL	1.50%	0.75	0.80%	0.4

TECHNICAL SPECIFICATION PIPE-TF-FIXED VOLUME PIPETTE						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
5µL	-	5µL	1.3%	0.065	1.2%	0.06
10µL	-	10µL	0.8%	0.08	0.8%	0.08
20µL	-	20µL	0.6%	0.12	0.5%	0.1
25µL	-	25µL	0.5%	0.125	0.3%	0.075
50µL	-	50µL	0.5%	0.25	0.3%	0.15
100µL	-	100µL	0.5%	0.5	0.3%	0.3
200µL	-	200µL	0.4%	0.8	0.2%	0.4
250µL	-	250µL	0.4%	1.0	0.2%	0.5
500µL	-	500µL	0.3%	1.5	0.2%	1.0
1000µL	-	1000µL	0.3%	3.0	0.2%	2.0
2000µL	-	2000µL	0.3%	6.0	0.15%	3.0
5000µL	-	5000µL	0.3%	15	0.15%	7.5



PIPE-MP, PIPE-MP8, PIPE-MP12, PIPE-MPF

FEATURES (SINGLE CHANNEL & MULTI-CHANNEL) DESCRIPTION

ITEM	DESCRIPTION
1	Fully autoclavable.
2	Ergonomic design provides excellent operating experience.
3	Easy-to-read volume display.
4	The pipettes cover the volume range from 0.1µL to 5mL.
5	Easy calibration and maintenance.
6	Manufactured from innovative materials.
7	Online calibration is available.

FEATURES (MULTI-CHANNEL) DESCRIPTION

ITEM	DESCRIPTION
1	Dispensing head rotates for effortless pipetting convenience.
2	Individual piston and tip cone assembly allow for easy repair and maintenance.
3	Compound material-made tip cone secures high sealing performance.
4	Compatible with most universal tip brands.
5	Online calibration is available.



PIPE - MP



PIPE - MP8



PIPE - MP12



PIPE - MPF

All rights reserved by iSiinNitzan. Reproduction and distribution of any part of this document is prohibited.

iSii-Nitzan GmbH

Via Cremignone 4C, Sorengo,
Switzerland, CH 6924

+(41) 765851837

info@isii-nitzan.swiss

www.isii-nitzan.swiss

More Information:



Versión: J2024-1

PIPE-MP, PIPE-MP8, PIPE-MP12, PIPE-MPF

TECHNICAL SPECIFICATION PIPE-MP SINGLE-CHANNEL ADJUSTABLE VOLUME PIPETTES						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
0.1-2.5µL	0.05µL	2.5µL	2.50%	0.0625	2.00%	0.05
		1.25µL	3.00%	0.0375	3.00%	0.0375
		0.25µL	12.00%	0.03	6.00%	0.015
0.5-10µL	0.1µL	10µL	1.00%	0.1	0.80%	0.08
		5µL	1.50%	0.075	1.50%	0.075
		1µL	2.50%	0.025	1.50%	0.015
2-20µL	0.5µL	20µL	0.90%	0.18	0.40%	0.08
		10µL	1.20%	0.12	1.00%	0.1
		2µL	3.00%	0.06	2.00%	0.04
5-50µL	0.5µL	50µL	0.60%	0.3	0.30%	0.15
		25µL	0.90%	0.225	0.60%	0.15
		5µL	2.00%	0.1	2.00%	0.1
10-100µL	1µL	100µL	0.80%	0.8	0.15%	0.15
		50µL	1.00%	0.5	0.40%	0.2
		10µL	3.00%	0.3	1.50%	0.15
20-200µL	1µL	200µL	0.60%	1.2	0.15%	0.3
		100µL	0.80%	0.8	0.30%	0.3
		20µL	3.00%	0.6	1.00%	0.2
50-200µL	1µL	200µL	0.60%	1.2	0.15%	0.3
		100µL	0.80%	0.8	0.30%	0.3
		50µL	1.00%	0.5	0.40%	0.2
100-1000µL	5µL	1000µL	0.60%	6	0.20%	2
		500µL	0.70%	3.5	0.25%	1.25
		100µL	2.00%	2	0.70%	0.7
200-1000µL	5µL	1000µL	0.60%	6	0.20%	2
		500µL	0.70%	3.5	0.25%	1.25
		200µL	0.90%	1.8	0.30%	0.6
1000-5000µL	50µL	5000µL	0.50%	25	0.15%	7.5
		2500µL	0.60%	15	0.30%	7.5
		1000µL	0.70%	7	0.30%	3
2-10mL	0.1mL	10mL	0.60%	60	0.20%	20
		5mL	1.20%	60	0.30%	15
		2mL	3.00%	60	0.60%	12

TECHNICAL SPECIFICATION PIPE-MP8 CHANNEL ADJUSTABLE VOLUME PIPETTES						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
0.5-10µL	0.1µL	10µL	1.50%	0.15	1.50%	0.15
		5µL	2.50%	0.125	2.50%	0.125
		1µL	4.00%	0.04	4.00%	0.04
5-50µL	0.5µL	50µL	1.00%	0.5	0.50%	0.25
		25µL	1.50%	0.375	1.00%	0.25
		5µL	3.00%	0.15	2.00%	0.1
50-300µL	5µL	300µL	0.70%	2.1	0.25%	0.75
		150µL	1.00%	1.5	0.50%	0.75
		50µL	1.50%	0.75	0.80%	0.4

TECHNICAL SPECIFICATION PIPE-MP12 CHANNEL ADJUSTABLE VOLUME PIPETTES						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
0.5-10µL	0.1µL	10µL	1.50%	0.15	1.50%	0.15
		5µL	2.50%	0.125	2.50%	0.125
		1µL	4.00%	0.04	4.00%	0.04
5-50µL	0.5µL	50µL	1.00%	0.5	0.50%	0.25
		25µL	1.50%	0.375	1.00%	0.25
		5µL	3.00%	0.15	2.00%	0.1
50-300µL	5µL	300µL	0.70%	2.1	0.25%	0.75
		150µL	1.00%	1.5	0.50%	0.75
		50µL	1.50%	0.75	0.80%	0.4

TECHNICAL SPECIFICATION PIPE-MPF FIXED VOLUME PIPETTE						
VOLUME RANGE	INCREMENT	TEST VOLUME	(ACCURACY ERROR)		(PRECISION ERROR)	
			%	µL	%	µL
5µL	-	5µL	1.3%	0.065	1.2%	0.06
10µL	-	10µL	0.8%	0.08	0.8%	0.08
20µL	-	20µL	0.6%	0.12	0.5%	0.1
25µL	-	25µL	0.5%	0.125	0.3%	0.075
50µL	-	50µL	0.5%	0.25	0.3%	0.15
100µL	-	100µL	0.5%	0.5	0.3%	0.3
200µL	-	200µL	0.4%	0.8	0.2%	0.4
250µL	-	250µL	0.4%	1.0	0.2%	0.5
500µL	-	500µL	0.3%	1.5	0.2%	1.0
1000µL	-	1000µL	0.3%	3.0	0.2%	2.0
2000µL	-	2000µL	0.3%	6.0	0.15%	3.0
5000µL	-	5000µL	0.3%	15	0.15%	7.5

All rights reserved by iSiinitzan. Reproduction and distribution of any part of this document is prohibited.

